


Welch Allyn and Cerner Vitals Connectivity



When it comes to vitals documentation, every second counts. That's why Welch Allyn and Cerner have partnered to help clinicians like you work more efficiently without compromising care.

Welch Allyn Connex® vital signs devices connect with the Cerner VitalsLink solution to send verified* vital signs data from the device at the bedside to the EMR—**no manual intervention required.**

A connected solution can help you improve clinical workflows and reduce errors.



Unconnected workflow:

- Take patient's vital signs
- Document on paper
- Find a workstation
- Open the patient's chart on a workstation and type vital signs into the EMR



Connected workflow:


- Scan the clinician ID badge
- Scan the patient ID
- Take patient's vital signs
- Tap "save." Vital signs are automatically sent to the EMR.

* Verification: Clinician scans patient wristband and reviews vitals before pressing save to send the results to the patient record.

Connectivity Counts

The impact on time savings, accuracy and clinical workflows is significant.

- **Save time** previously spent manually printing, scanning and importing data into the EMR
- **Improve accuracy** by avoiding manual processes and associated errors
- **Work more efficiently** by freeing up time for other priorities



The following Welch Allyn vital signs devices all interface with Cerner: Connex Spot Monitor, Connex Vital Signs Monitor and Connex Integrated Wall System. Contact your Welch Allyn representative for more information on how to connect.